




Radiographic Fluoroscopy C Arm DRF Machine 380V With High Frequency Generator

Our Product Introduction

Basic Information

- Place of Origin: CHINA
- Brand Name: VOLITION
- Certification: ISO13485、FSC
- Model Number: Picke 2800DRF
- Minimum Order Quantity: 1
- Price: Negotiable
- Packaging Details: Standard Exportation Package
- Delivery Time: 7-15 work days
- Payment Terms: LC, T/T, D/P, PayPal, Western Union




Product Specification

- Type: Floor-mounted
- Power Supply: 110V/220V/380V
- Output Power: 50kw High Frequency Generator
- X-ray Tube: Cannon
- Column: Motor-driven
- Detector Size: 17"x 17"
- Detector Type: Dynamic
- Highlight: Radiographic Fluoroscopy DRF Machine, 380V DRF Machine, 380V c arm fluoroscopy machine



More Images



for more products please visit us on xray-system.com

Product Description

The patented multifunctional lifting bed can perform retrograde urography; Hysterosalpingography; T-tube imaging; Percutaneous

cholangiography;

Bedside instrument tray for convenient intervention operations.

Bedside proximity control panel: capable of intelligent control of various movements; Generator parameter adjustment; And can perform image

control.

Bedside display: convenient for nearby operators to observe images in real-time

Rotation angle range: $+90^{\circ} \sim -15^{\circ}$; Bed lifting range: $\geq 400\text{mm}$;

Horizontal movement range of bed surface: $100\text{mm} \pm 10\%$;

X-ray tube lifting range: $\geq 1200\text{mm}$; Ball tube rotation range: $\pm 90^{\circ}$;

The bed rotates, the bed rises and falls, the bed surface moves, and the detector moves within its range of motion, allowing for individual or

composite movements at will; Electromagnetic locking.

*Customizable according to customer requirements



Liuyang Volition Medical Equipment Co., Ltd.



+86 17775818268



258031353@qq.com



xray-system.com

76 Jianshou Ave, Economic and Technological Development Zone, Liuyang, Hunan, China