

## OEM Mobile C Arm Xray Dental Radiology Equipment 50kw

Our Product Introduction

### Basic Information

- Place of Origin: CHINA
- Brand Name: VOLITION
- Certification: ISO13485、FSC
- Model Number: Picke 2200DRF
- Minimum Order Quantity: 1
- Price: Negotiable
- Packaging Details: Standard Exportation Package
- Delivery Time: 7-15 work days
- Payment Terms: LC, T/T, D/P, PayPal, Western Union



### Product Specification

- Application: Metal Processing
- Power Supply: 110V/220V
- Type: Mobile C-arm
- X-ray Tube Heat Capacity: 300Khu (210KJ)
- Dimensions: 1000mm X 500mm X 800mm
- Max Photograph Tube Current: 640mA
- Warranty: 2 Years
- Motorized C-arm: Yes
- Max Output Power: 50kw
- Highlight: **Mobile dental radiology equipment,  
50kw dental radiology equipment,  
OEM mobile xray machine**

for more products please visit us on [xray-system.com](http://xray-system.com)

## Product Description

Excellent image algorithm without dynamic tailing

AK & AKR dynamic monitoring; Real time dose display; High dose prompt; Over dose suspension exposure;

50KW high power.

Simple structure, highly integrated, no external drag cable, easy to use;

Large anode heat capacity to support long-term use



**Liuyang Volition Medical Equipment Co., Ltd.**



+86 17775818268



258031353@qq.com



xray-system.com

76 Jianshou Ave, Economic and Technological Development Zone, Liuyang, Hunan, China