




50kw Digital Fluoroscopy Machine Radiography Medical X Ray Equipment

Our Product Introduction

Basic Information

- Place of Origin: CHINA
- Brand Name: VOLITION
- Certification: ISO13485、FSC
- Model Number: Picke 1800DR/500
- Minimum Order Quantity: 1
- Price: Negotiable
- Packaging Details: Standard Exportation Package
- Delivery Time: 7-15 work days
- Payment Terms: LC, T/T, D/P, PayPal, Western Union



Product Specification

- Max Tube Current: 500mA
- Detector Pixel Matrix: 3072 X 3072
- Tube Heat Capacity: 200KHU
- Diagnostic Display: Monochrome
- Power Supply: 110V/220V
- Exposure Time: 0.001 - 6 Seconds
- Max Output Power: 50kw
- Current Time Product: 0.1mAs~630mAs
- Detector And Tube Automatic Alignment Function: Yes
- Highlight: 50kw digital fluoroscopy machine,
XRay digital fluoroscopy machine,
Fluoroscopy medical x ray equipment

for more products please visit us on xray-system.com

Product Description

50Kw high-frequency generator and high heat capacity X-ray tube provide strong guarantee for outputting high-quality X-rays;

Adopting the company's original image acquisition system and a powerful image algorithm independently developed by the R&D team; Adhering to

the principle of image first, ensuring clear images;

Electric chest X-ray holder with automatic tracking and alignment and compartment remote control function, easy and simple to operate; Can

connect to hospital PACS system; Remote maintenance and parameter calibration; The image post-processing function is powerful.



Liuyang Volition Medical Equipment Co., Ltd.



+86 17775818268



258031353@qq.com



xray-system.com

76 Jianshou Ave, Economic and Technological Development Zone, Liuyang, Hunan, China