



## FSC Digital Mobile Medical Radiology Equipment Fluoroscope X Ray Machines 50kw

Our Product Introduction

### Basic Information

- Place of Origin: CHINA
- Brand Name: VOLITION
- Certification: ISO13485、FSC
- Model Number: PICKE 1900 DRM
- Minimum Order Quantity: 1
- Price: Negotiable
- Packaging Details: Standard Exportation Package
- Delivery Time: 7-15 work days
- Payment Terms: LC, T/T, D/P, PayPal, Western Union



### Product Specification

- Anode Heat Capacity: 150KHU
- Max Output Power: 50kw
- Material: Stainless Steel
- Image Resolution:  $\geq 3.5\text{lp/mm}$
- Products: Mobile Medical X Ray Machine
- Max Current: 500mA
- Tube Voltage: 40-150KV
- Loading Time: 0.001s-6.3s
- FPD: 17"\*17" Wireless
- Screen: Touch Screen
- Highlight: Fluoroscope medical radiology equipment,  
FSC medical radiology equipment,  
FSC fluoroscope x ray machines

for more products please visit us on [xray-system.com](http://xray-system.com)

## Product Description

Pick-1900 DRM is an economical and practical digital X-ray system, which can be widely used in radiology, orthopedics, health examination,

emergency room, ICU, operating room and neonatology department. Fully meet the needs of digital radiation in hospital.

40kW super power high frequency and high voltage generator and high heat capacity ball tube provide powerful guarantee for high quality X-ray

output. Wireless flat panel detector, get rid of cable bondage; Wide nose angle, flexible positioning, hovering at any position within the angle range; The battery has a large energy storage capacity.



**Liuyang Volition Medical Equipment Co., Ltd.**



+86 17775818268



258031353@qq.com



xray-system.com

76 Jianshou Ave, Economic and Technological Development Zone, Liuyang, Hunan, China