

Mobile DR Digital Radiography Fluoro X Ray Machine 50kw

Our Product Introduction

for more products please visit us on xray-system.com

Basic Information

- Place of Origin: CHINA
- Brand Name: VOLITION
- Certification: ISO13485、FSC
- Model Number: PICKE 1900 DRM
- Minimum Order Quantity: 1
- Price: Negotiable
- Packaging Details: Standard Exportation Package
- Delivery Time: 7-15 work days
- Payment Terms: LC, T/T, D/P, PayPal, Western Union



Product Specification

- Features: Easy To Maneuver And Positioning
- Dimensions: Compact
- Type: Mobile X Ray Machine
- Display: Touchscreen LCD
- Image Quality: High Resolution
- Technology: Digital Radiography
- Detector Type: Wireless FPD
- Power Supply: 110V/220V
- Battery: Yes
- Highlight: **DR fluoro x ray machine,
50kw fluoro x ray machine,
DR digital radiography x ray machine**



Product Description

Pick-1900 DRM is an economical and practical digital X-ray system, which can be widely used in radiology, orthopedics, health examination,

emergency room, ICU, operating room and neonatology department. Fully meet the needs of digital radiation in hospital.

40kW super power high frequency and high voltage generator and high heat capacity ball tube provide powerful guarantee for high quality X-ray

output. Wireless flat panel detector, get rid of cable bondage; Wide nose angle, flexible positioning, hovering at any position within the angle range; The battery has a large energy storage capacity.



Liuyang Volition Medical Equipment Co., Ltd.



+86 17775818268



258031353@qq.com



xray-system.com

76 Jianshou Ave, Economic and Technological Development Zone, Liuyang, Hunan, China