



## Lightweight Portable X Ray Equipment Veterinary Radiology Equipment 220V

Our Product Introduction

### Basic Information

- Place of Origin: China
- Brand Name: VOLITION
- Certification: ISO13485, FSC
- Minimum Order Quantity: 1
- Price: Negotiable
- Packaging Details: Standard Exportation Package
- Delivery Time: 7-15 work days
- Payment Terms: LC, T/T, D/P, PayPal, Western Union



### Product Specification

- Image Receptor Type: Flat Panel Detector
- Weight: Less Than 10 Kg
- Usage: Medical Imaging
- Focal Spot Size: 1.5 Mm
- Image Receptor Size: 14x17 Inches
- Exposure Time: 0.01-5 Seconds
- Image Resolution: 3.4 Lp/mm
- Power Source: Battery
- Highlight: Radiology portable x ray equipment,  
220V portable x ray equipment,  
Portable veterinary radiology equipment

for more products please visit us on [xray-system.com](http://xray-system.com)

## Product Description

### Product Description:

Whether you are looking for a 3D dental x ray machine or digital radiography equipment for other medical imaging applications, this portable radiography machine is an excellent choice. With its large 17"x17" detector size and advanced imaging capabilities, it is ideal for a wide range of diagnostic imaging tasks, including dental x ray equipment and more.

One of the key benefits of this portable radiography machine is its lightweight and compact design, which makes it easy to transport and use in a variety of settings. Whether you are working in a hospital, clinic, or other medical facility, this machine can be quickly set up and used to capture high-quality x-ray images with ease.

Other features of this portable radiography machine include a range of imaging modes and settings, as well as advanced image processing capabilities that help to enhance the quality and clarity of your x-ray images. Whether you are using this machine for routine diagnostic imaging tasks or more complex medical imaging procedures, it is designed to provide the performance and flexibility you need. In conclusion, the Portable Radiography Machine is an exceptional medical imaging device that offers advanced features and capabilities in a portable and easy-to-use package. Whether you are looking for a 3D dental x ray machine, dental x ray equipment, or digital radiography equipment for other medical imaging applications, this machine is an excellent choice. With its large image receptor size, high X-ray tube current, and advanced imaging capabilities, it is designed to provide exceptional performance and flexibility for a wide range of diagnostic imaging tasks.

### Features:

Product Name: Portable Radiography Machine

Detector size: 17"x17"

X-Ray Tube Current: 5 MA

Focal Spot Size: 1.5 Mm

Power Supply: 110V/220V/380V

Usage: Medical Imaging

This Portable Radiography Machine is a state-of-the-art digital radiography equipment that can be used for both veterinary and dental purposes. With its 17"x17" detector size, it can capture high-quality images with a 5 MA X-Ray Tube Current and 1.5 mm focal spot size. This portable radiography machine can be powered by 110V, 220V, or 380V, making it suitable for use in various settings. It is a great alternative to conventional radiography machines and is a must-have for any medical imaging practice.

### Technical Parameters:

Power Source:	Battery
Usage:	Medical Imaging
X-Ray Tube Type:	Stationary
Weight:	Less Than 10 Kg
Image Receptor Size:	14x17 Inches
Image Receptor Type:	Flat Panel Detector
X-Ray Tube Current:	5 MA
Collimator:	LED
Image Resolution:	3.4 Lp/mm
Type:	Portable

### Applications:

The VOLITION Portable Radiography Machine is equipped with a 40kw/32kw high-frequency generator that provides high-quality images with minimal radiation exposure. The machine's focal spot size is 1.5 mm, which is ideal for capturing detailed images of different body parts. It also has a flat panel detector as an image receptor type, which is known for its high resolution and accuracy.

The VOLITION Portable Radiography Machine is suitable for use in various medical scenarios. For instance, it can be used in emergency departments, intensive care units, and operating rooms. It is also an excellent choice for dental clinics, where it can be used as a 3d dental x ray machine to capture detailed images of teeth and gums.

The VOLITION Portable Radiography Machine is designed to improve workflow efficiency and reduce patient waiting times. Its portability allows it to be easily moved to the patient's location, eliminating the need for patients to be transported to a separate imaging room. This feature is particularly useful in emergency situations where time is of the essence.

In conclusion, the VOLITION Portable Radiography Machine is a versatile and high-quality machine that is ideal for medical imaging in various settings. Its lightweight and portable design, combined with its high-quality imaging capabilities, make it an excellent choice for healthcare providers looking for a reliable and efficient radiography machine.

### Customization:

### Packing and Shipping:

#### Product Packaging:

Portable Radiography Machine

Power cord

Remote control

User manual

Protective carrying case

#### Shipping:

Shipping method: UPS ground

Expected delivery time: 3-5 business days

Shipping cost: \$50

Tracking number will be provided



**Liuyang Volition Medical Equipment Co., Ltd.**



+86 17775818268



258031353@qq.com



xray-system.com

76 Jianshou Ave, Economic and Technological Development Zone, Liuyang, Hunan, China