

OEM Digital Radiography X Ray Machine System For Dental Medical Imaging

Our Product Introduction

for more products please visit us on xray-system.com

Basic Information

- Place of Origin: China
- Brand Name: VOLITION
- Certification: ISO13485、FSC
- Minimum Order Quantity: 1
- Price: Negotiable
- Packaging Details: Standard Exportation Package
- Delivery Time: 7-15 work days
- Payment Terms: LC, T/T, D/P, PayPal, Western Union



Product Specification

- X-Ray Tube Current: 5 MA
- Power Supply: 110V/220V/380V
- Exposure Time: 0.01-5 Seconds
- Weight: Less Than 10 Kg
- Detector Size: 17"x17"
- Type: Portable
- Image Receptor Size: 14x17 Inches
- Usage: Medical Imaging
- Highlight: **OEM digital radiography x ray machine,
Dental digital radiography x ray machine,
OEM digital x ray system**

Product Description

Product Description:

The Portable Radiography Machine is equipped with a powerful battery power source, allowing users to use the machine without any external power source. This feature makes it an ideal choice for emergency situations where power sources may not be readily available. With its lightweight design, the Portable Radiography Machine weighs less than 10 kg, making it easy to carry around and use in different locations.

The Portable Radiography Machine features a 17"x17" detector size, which is perfect for capturing detailed images of different body parts. This feature ensures that medical professionals can make accurate diagnoses and provide appropriate treatments to patients. The Portable Radiography Machine is suitable for a wide range of medical imaging applications, including orthopedics, gastroenterology, cardiology, and more.

The Portable Radiography Machine is easy to operate and features a user-friendly interface. The machine comes with intuitive controls that allow users to quickly and easily capture images. The Portable Radiography Machine is designed to be durable and reliable, ensuring that it can withstand daily use in demanding medical environments.

The Portable Radiography Machine is a versatile and efficient medical imaging device that is perfect for medical professionals who require mobility and convenience. This machine is an ideal choice for emergency medical situations, where quick and accurate diagnoses are critical. The Portable Radiography Machine is a top-of-the-line handheld x ray machine that offers excellent performance and reliability, making it an excellent investment for any medical facility.

Features:

Product Name: Portable Radiography Machine

X-ray Tube: Cannon

Image Resolution: 3.4 Lp/mm

Exposure Time: 0.01-5 Seconds

Power Supply: 110V/220V/380V

Image Receptor Type: Flat Panel Detector

This new portable x-ray machine is perfect for industrial radiography equipment needs. It features a Cannon x-ray tube and an image resolution of 3.4 Lp/mm. The exposure time ranges from 0.01-5 seconds and it can be powered by 110V, 220V, or 380V. The image receptor type is a flat panel detector, making it a modern and computerised x-ray machine.

Technical Parameters:

Power Supply	110V/220V/380V
Weight	Less Than 10 Kg
Detector Size	17"x17"
Collimator	LED
Type	Portable
Usage	Medical Imaging
Focal Spot Size	1.5 Mm
Power Source	Battery
Image Resolution	3.4 Lp/mm
Image Receptor Type	Flat Panel Detector

Applications:

One of the key features of the VOLITION Portable Radiography Machine is its high image resolution of 3.4 Lp/mm. This ensures that the images produced are clear and detailed, making it an ideal choice for a wide range of radiology equipment applications. The focal spot size of 1.5 mm also ensures that the images produced are sharp and precise, making it an excellent choice for critical diagnostic procedures.

Despite its impressive capabilities, the VOLITION Portable Radiography Machine is also incredibly lightweight, weighing less than 10 kg. This makes it easy to transport and use in a variety of different settings, from hospitals and clinics to remote or mobile locations.

The image receptor size of 14x17 inches is another key feature of this portable x-ray machine. This size is suitable for a wide range of imaging applications, including chest X-rays, abdominal X-rays, and more. The machine's portable design also means that it can be used in a variety of different scenarios, from emergency medical situations to veterinary clinics and more.

Overall, the VOLITION Portable Radiography Machine is an excellent choice for anyone looking for a high-quality, portable x-ray machine that can be used in a variety of different settings. Whether you are a healthcare professional looking for a reliable radiology equipment option, or a veterinarian in need of a mobile x-ray machine, the VOLITION Portable Radiography Machine is an excellent choice.

Customization:

Customize your VOLITION portable radiography machine to fit your specific needs:

Brand Name: VOLITION

Place of Origin: China

Certification: ISO13485,FSC

Power Source: Battery

Collimator: LED

Image Resolution: 3.4 Lp/mm

Image Receptor Type: Flat Panel Detector

X-Ray Tube Voltage: 70 KV

Our customization services ensure that your new portable x ray machine meets all of your dental x ray machine needs. With digital radiography equipment capabilities, you can trust your VOLITION portable radiography machine to deliver high-quality images every time.

Packing and Shipping:

Product Packaging:

The Portable Radiography Machine will be carefully packaged in a sturdy cardboard box with ample padding to prevent any damage during shipping. The machine will be secured with foam inserts and wrapped in plastic to protect it from any moisture or dust. The box will be clearly labeled with the product name, model number, and handling instructions.

Shipping:

The Portable Radiography Machine will be shipped via a reputable courier service, such as FedEx or UPS, to ensure timely and safe delivery. The shipping cost will be calculated based on the destination and weight of the package. Customers will receive a tracking number so they can follow their package's progress and know when to expect delivery.



Liuyang Volition Medical Equipment Co., Ltd.



+86 17775818268



258031353@qq.com



xray-system.com

76 Jianshou Ave, Economic and Technological Development Zone, Liuyang, Hunan, China