


ISO13485 Digital Portable Radiography Machine Veterinary X Ray System

Our Product Introduction

for more products please visit us on xray-system.com

Basic Information

- Place of Origin: China
- Brand Name: VOLITION
- Certification: ISO13485、FSC
- Minimum Order Quantity: 1
- Price: Negotiable
- Packaging Details: Standard Exportation Package
- Delivery Time: 7-15 work days
- Payment Terms: LC, T/T, D/P, PayPal, Western Union



Product Specification

- Detector Size: 17"x17"
- X-Ray Tube Type: Stationary
- Power Source: Battery
- X-Ray Tube Voltage: 70 KV
- Exposure Time: 0.01-5 Seconds
- Type: Portable
- Focal Spot Size: 1.5 Mm
- Image Receptor Type: Flat Panel Detector
- Highlight: **ISO13485 Portable Radiography Machine, Veterinary XRay Portable Radiography Machine, ISO13485 portable digital x ray machine**

Product Description

Product Description:

One of the most remarkable features of this Portable Radiography Machine is its versatility in use. It can be used for dental x ray equipment, as well as 3d dental x ray machine, making it a highly valuable device for dentists and orthodontists. Its compact size and lightweight design make it easy to move around the clinic or hospital, and it can be easily stored when not in use.

The high frequency generator in this Portable Radiography Machine produces consistent and reliable power output, which results in clear and accurate images. It also reduces the exposure time for patients, making the procedure faster and more comfortable. The machine is also equipped with advanced safety features to ensure that both patients and technicians are protected from radiation exposure.

This Portable Radiography Machine is an excellent investment for any medical facility looking to upgrade their medical imaging capabilities. It is ideal for hospitals, clinics, and mobile healthcare units, and it provides an excellent return on investment. The machine is durable and easy to maintain, ensuring that it will last for many years.

In conclusion, if you are looking for a versatile, reliable, and efficient Portable Radiography Machine for your medical facility, then this product is an excellent choice. With its advanced features, it is perfect for dental x ray equipment, 3d dental x ray machine, and other medical imaging applications. Get yours today and experience the benefits of high-quality digital x-ray system.

Features:

Product Name: Portable Radiography Machine

Image Receptor Size: 14x17 Inches

Image Resolution: 3.4 Lp/mm

Detector size: 17"x17"

Output Power: 40kw/32kw High Frequency Generator

X-Ray Tube Type: Stationary

Can be used as a handheld x ray machine

Mobile xray machine, easy to transport

Suitable for veterinary dental x ray machine purposes

Technical Parameters:

Power Source	Battery
Focal Spot Size	1.5 Mm
Image Resolution	3.4 Lp/mm
Exposure Time	0.01-5 Seconds
Collimator	LED
X-Ray Tube Voltage	70 KV
X-ray Tube	Cannon
Power Supply	110V/220V/380V
Detector size	17"x17"
X-Ray Tube Current	5 MA

This technical parameter table is for a portable radiography machine, which can be used as a dental x ray machine or industrial radiography equipment. It has a battery power source, LED collimator, and a Cannon x-ray tube. The machine has a 70 KV x-ray tube voltage, 5 MA x-ray tube current, and 17"x17" detector size. The image resolution is 3.4 Lp/mm and the exposure time can range from 0.01-5 seconds. The power supply options are 110V/220V/380V.

Applications:

Thanks to its powerful 40kw/32kw High Frequency Generator, this digital radiography equipment is capable of producing crystal-clear images that are easy to interpret and diagnose. Its large 17"x17" detector size and 14x17 inch Image Receptor Size make it ideal for a wide range of medical applications, from orthopedic surgery to cardiac imaging.

One of the biggest advantages of the VOLITION Portable Radiography Machine is its versatility. Whether you're working in a hospital, clinic, or private practice, this digital x ray machine is easy to transport and set up, allowing you to take high-quality images wherever and whenever you need them. Its battery power source also ensures that you'll never be tethered to an electrical outlet.

Some of the most common application occasions and scenarios for the VOLITION Portable Radiography Machine include:

Emergency rooms and trauma centers, where fast and accurate imaging is essential for saving lives.

Operating rooms, where surgeons need to be able to quickly and easily visualize the area they're operating on.

ICUs and CCUs, where patients with serious medical conditions require constant monitoring and imaging.

Private practices and clinics, where patients can benefit from the convenience and flexibility of on-site imaging services.

Overall, the VOLITION Portable Radiography Machine is a powerful and reliable piece of equipment that is well-suited for a wide range of medical applications. Whether you're a doctor, nurse, or other medical professional, this digital radiography equipment can help you provide better care to your patients and make more accurate diagnoses.

Customization:

Packing and Shipping:

Product Packaging:

The Portable Radiography Machine will be securely packed and shipped in a sturdy cardboard box with appropriate cushioning to prevent any damage during transport. The box will also be labeled with handling instructions and marked as fragile to ensure proper handling.

Shipping:

The Portable Radiography Machine will be shipped via a reputable courier service, and the shipping cost will be included in the purchase price. The estimated delivery time will depend on the destination and the courier's schedule, but we will provide tracking information so you can monitor the progress of your shipment.



Liuyang Volition Medical Equipment Co., Ltd.



+86 17775818268



258031353@qq.com



xray-system.com

76 Jianshou Ave, Economic and Technological Development Zone, Liuyang, Hunan, China