

C Arm Dental Xray Unit Digital Radiography Equipment FSC Certificated

Our Product Introduction

for more products please visit us on xray-system.com

Basic Information

- Place of Origin: China
- Brand Name: VOLITION
- Certification: ISO13485、FSC
- Minimum Order Quantity: 1
- Price: Negotiable
- Packaging Details: Standard Exportation Package
- Delivery Time: 7-15 work days
- Payment Terms: LC, T/T, D/P, PayPal, Western Union



Product Specification

- Focal Spot Size: 1.5 Mm
- Weight: Less Than 10 Kg
- Detector Size: 17"x17"
- Exposure Time: 0.01-5 Seconds
- Power Source: Battery
- Image Resolution: 3.4 Lp/mm
- X-Ray Tube Type: Stationary
- X-ray Tube: Cannon
- Highlight: **C arm dental xray unit,
FSC certificated dental xray unit,
Dental digital radiography equipment**

Product Description

Product Description:

Equipped with a powerful X-Ray Tube Current of 5 MA, the Portable Radiography Machine delivers clear and accurate images that allow clinicians to make informed decisions about patient care. The machine also features a LED collimator that enables precise positioning and reduces unnecessary radiation exposure.

Powered by a battery, the Portable Radiography Machine can be used anywhere, even in remote locations without access to electrical outlets. This makes it an ideal choice for medical professionals who need to take x-rays in the field or in emergency situations.

The Portable Radiography Machine utilizes a Flat Panel Detector as its image receptor type, which provides fast and efficient image capture with minimal radiation exposure. The images produced by this machine are of the highest quality and can be easily shared with other medical professionals for consultation or review.

In summary, the Portable Radiography Machine is a state-of-the-art handheld x ray machine that delivers powerful diagnostic capabilities in a portable, easy-to-use package. Whether you're working in a hospital, clinic, or out in the field, this digital x ray machine is the perfect tool for capturing high-quality images quickly and accurately.

Features:

Product Name: New Portable X Ray Machine

Image Receptor Size: 14x17 Inches

Power Source: Battery

Collimator: LED

Weight: Less Than 10 Kg

Usage: Medical Imaging

Features:

Computerised Xray Machine

Digital X Ray Machine

Technical Parameters:

X-Ray Tube Current	5 MA
Image Resolution	3.4 Lp/mm
Weight	Less Than 10 Kg
Collimator	LED
X-Ray Tube Voltage	70 KV
Type	Portable
Usage	Medical Imaging
Image Receptor Size	14x17 Inches
Image Receptor Type	Flat Panel Detector
Focal Spot Size	1.5 Mm

This is a table of technical parameters for a mobile xray machine, also known as a portable x ray machine or digital radiography equipment.

Applications:

C-arm Radiology Equipment: This portable x-ray unit is ideal for use in C-arm radiology equipment, allowing medical professionals to get clear, high-quality images of a patient's bones or organs during surgery. With an exposure time of 0.01-5 seconds and an X-ray tube current of 5 MA, this portable radiography machine delivers fast and accurate results.

Emergency Medical Services: In emergency medical services, speed and efficiency are crucial. The VOLITION portable radiography machine can be quickly transported to the patient's location, allowing medical professionals to get a clear view of the injury or condition and begin treatment immediately.

Remote Locations: In remote locations where medical facilities are not readily available, this portable x-ray machine can be a lifesaver. It can be used to diagnose and treat injuries and illnesses in areas where traditional medical equipment may not be accessible.

Nursing Homes: Nursing homes can also benefit from having a portable radiography machine on hand. It allows medical professionals to quickly and easily diagnose and treat conditions without having to transport patients to a hospital or clinic.

With a power supply of 110V/220V/380V, this portable radiography machine is versatile and can be used in a variety of settings. Whether it's for emergency medical services, nursing homes, remote locations, or c-arm radiology equipment, the VOLITION portable x-ray machine is a reliable and efficient solution.

Customization:

Customization Services:

Brand Name: VOLITION

Place of Origin: China

Certification: ISO13485,FSC

Image Receptor Size: 14x17 Inches

X-Ray Tube Voltage: 70 KV

X-Ray Tube Type: Stationary

Power Supply: 110V/220V/380V

Power Source: Battery

Our Portable Radiography Machine product can be customized to fit your specific needs. Whether you need a conventional radiography machine, mobile xray machine, or a handheld x ray machine, we can make it happen. Contact us for more information on our customization services.

Packing and Shipping:

Product Packaging:

The Portable Radiography Machine will be carefully packaged to ensure safe delivery. The machine will be wrapped in protective material and placed in a sturdy box with additional padding to prevent any damage during shipping.

Shipping:

Shipping will be handled by our trusted delivery partners. We offer free standard shipping for all orders within the United States. For international orders, additional shipping fees may apply. We will provide a tracking number once the order has shipped so you can easily track the progress of your delivery.



Liuyang Volition Medical Equipment Co., Ltd.



+86 17775818268



258031353@qq.com



xray-system.com

76 Jianshou Ave, Economic and Technological Development Zone, Liuyang, Hunan, China