

## Customized Veterinary Digital Radiography System X Ray Fluoroscopy Machine

Our Product Introduction

for more products please visit us on [xray-system.com](http://xray-system.com)

### Basic Information

- Place of Origin: China
- Brand Name: VOLITION
- Certification: ISO13485, FSC
- Minimum Order Quantity: 1
- Price: Negotiable
- Packaging Details: Standard Exportation Package
- Delivery Time: 7-15 work days
- Payment Terms: LC, T/T, D/P, PayPal, Western Union



### Product Specification

- X-Ray Generator: High Frequency Or Constant Potential
- X-Ray Tube Focal Spot: 0.3-1.5 Mm
- Type: Fixed Or Mobile
- Power Supply: AC Or DC
- Image Intensifier Field Of View: 9 Inch: 23 Cm, 12 Inch: 30 Cm
- Tube Voltage: 50-150 KV
- Tube Current: 0.5-5 MA
- Image Intensifier: 9 Inch Or 12 Inch
- Highlight:
  - customized veterinary digital radiography system
  - Fluoroscopy veterinary digital radiography system
  - Radiography x ray fluoroscopy machine

## Product Description

### Product Description:

Our Fluoroscope X Ray Machine also comes with a high-frequency or constant potential x-ray generator for optimal performance. With a tube current range of 0.5-5 MA, our machine offers versatility for all your imaging needs. Whether you're in the medical industry or a veterinary practice, our digital x ray machine veterinary is perfect for your needs. Our machine offers the highest quality images for accurate diagnoses. Our mobile digital x ray machine is also available for those on-the-go imaging needs. Our machine is portable and easy to move, making it perfect for field work or emergency situations. Trust in our Fluoroscope X Ray Machines for all your imaging needs. Contact us today to learn more about our product and how it can benefit your industry.

### Features:

Product Name: Fluoroscope X Ray Machines  
Image Intensifier Magnification: 1.5-2.5x  
X-Ray Generator: High Frequency Or Constant Potential  
Image Processing: Digital Or Analog  
Image Intensifier Material: Cesium Iodide Or Gadolinium Oxysulfide  
Image Intensifier Field Of View: 9 Inch: 23 Cm, 12 Inch: 30 Cm

### Technical Parameters:

Image Intensifier Field Of View	9 Inch: 23 Cm, 12 Inch: 30 Cm
Image Intensifier Resolution	1.5-2.5 Lp/mm
Tube Voltage	50-150 KV
X-Ray Tube Focal Spot	0.3-1.5 Mm
X-Ray Tube Cooling	Air Or Oil
Image Processing	Digital Or Analog
Image Intensifier Material	Cesium Iodide Or Gadolinium Oxysulfide
Image Intensifier	9 Inch Or 12 Inch
Tube Current	0.5-5 MA
X-Ray Generator	High Frequency Or Constant Potential

Our medical x ray equipment includes a digital x ray machine veterinary and veterinary digital dental radiography equipment.

### Applications:

Image Intensifier: 9 Inch Or 12 Inch  
X-Ray Tube Cooling: Air Or Oil  
Image Processing: Digital Or Analog  
Power Supply: AC Or DC  
Type: Fixed Or Mobile

The VOLITION Fluoroscope X Ray Machine is suitable for a wide range of medical applications, and some of the common product application occasions and scenarios include:

**Hospital Radiology Department:** The VOLITION Fluoroscope X Ray Machine is ideal for use in hospital radiology departments for a wide range of imaging applications. The machine can be configured for digital or analog image processing, and it is available in both fixed and mobile configurations to meet different requirements.

**Small Portable X Ray Machine:** The VOLITION Fluoroscope X Ray Machine is also available in a compact, portable configuration that is ideal for use in small clinics and mobile medical units. The small portable x ray machine can be easily transported from one location to another, and it is designed to provide high-quality imaging in a convenient and efficient manner.

**Veterinary Xray Equipment:** The VOLITION Fluoroscope X Ray Machine is also suitable for use in veterinary clinics and hospitals. The machine can be configured to meet the specific needs of veterinary applications, and it is available with a variety of image intensifier sizes and x-ray tube cooling options to ensure optimal performance.

**Portable Dental Xray Unit:** The VOLITION Fluoroscope X Ray Machine is also ideal for use in dental clinics and mobile dental units. The machine can be configured for digital or analog image processing, and it is available in both fixed and mobile configurations to meet different requirements. The compact and portable design of the machine makes it easy to transport and use in different locations.

Overall, the VOLITION Fluoroscope X Ray Machine is a versatile and reliable medical imaging equipment that can support a wide range of applications. Whether you need a fixed or mobile machine, digital or analog image processing, or a specific image intensifier size or x-ray tube cooling option, the VOLITION Fluoroscope X Ray Machine can be configured to meet your specific needs.

## Customization:

## Packing and Shipping:

### Product Packaging for Fluoroscope X Ray Machines:

The machine will be securely packaged in a wooden crate to prevent any damage during shipping.

All accessories, manuals, and other necessary documents will be included in the package.

The package will be labeled with all the necessary information, including the recipient's address, the weight, and the dimensions of the crate.

### Shipping:

The package will be shipped via a reputable courier service, such as FedEx or DHL.

The shipping method will be chosen based on the customer's preferences and the delivery timeframe.

The shipping cost will be calculated based on the weight and dimensions of the package, as well as the destination.

The customer will be provided with a tracking number to monitor the status of the shipment.



**Liuyang Volition Medical Equipment Co., Ltd.**



+86 17775818268



258031353@qq.com



xray-system.com

76 Jianshou Ave, Economic and Technological Development Zone, Liuyang, Hunan, China