

ODM Veterinary Dental Digital Fluoroscopy X Ray Machine High Precision Imaging

Our Product Introduction

for more products please visit us on xray-system.com

Basic Information

- Place of Origin: China
- Brand Name: VOLITION
- Certification: ISO13485, FSC
- Minimum Order Quantity: 1
- Price: Negotiable
- Packaging Details: Standard Exportation Package
- Delivery Time: 7-15 work days
- Payment Terms: LC, T/T, D/P, PayPal, Western Union



Product Specification

- Resistance: Not Larger Than 1Ω
- Kv Range: 40-150 KV
- Intermittent Mode: 100kV ,200mA,20kW
- Power Supply: 220V/50Hz
- Material: Stainless Steel
- Item: Portable X Ray Machine
- Rectifier Method: Full-wave Rectifier
- Washing Control: To Enter Water 2~8L/min Automatically
- Highlight: **ODM fluoroscopy x ray machine ,
Veterinary fluoroscopy x ray machine,
ODM digital fluoroscopy machine**

Product Description

Product Description:

One of the key features of this x ray machine is its KV range, which is 40-150 KV. This range allows for a wide range of diagnostic applications, making it suitable for use in a variety of medical settings. Additionally, this x ray machine has a resistance of not larger than 1Ω, ensuring that it is safe and reliable to use.

The Medical X Ray Machine also has a maximum current of 500mA, allowing for high-quality diagnostic images with minimal exposure time. This makes it a suitable option for busy medical facilities that need to perform a high volume of diagnostic exams in a short period of time.

Another important feature of this x ray machine is its intermittent mode, which is 100kV, 200mA, and 20kW. This mode allows for a more efficient use of energy and reduces the risk of overheating, making it a safer and more reliable option for medical professionals.

The Medical X Ray Machine is made of high-quality stainless steel, making it durable and long-lasting. This material is also easy to clean and maintain, ensuring that the x ray machine remains in good condition for years to come.

This x ray machine is suitable for a variety of medical applications, including handheld dental x ray machine, veterinary dental x ray machine, and other diagnostic imaging procedures.

In summary, the Medical X Ray Machine is a reliable and efficient tool for medical professionals in need of high-quality diagnostic images. Its wide KV range, maximum current, and intermittent mode make it a suitable option for busy medical facilities, while its durable stainless steel construction ensures that it will last for years to come.

Features:

Product Name: Medical X Ray Machine

Max: 0.4-360mAs

Kv Range: 40-150 KV

Intermittent Mode: 100kV, 200mA, 20kW

Nominal Power: 100kV, 40mA, 0.1s, 4kW

Max Current: 500mA

Can be used as dental xray unit

Can be used as veterinary radiology equipment

Equipped with x ray machine digital radiography

Technical Parameters:

Item	Portable X Ray Machine
Material	Stainless Steel
X Ray Tube Type	Rotating Anode
Voltage Frequency	380V/220V 50HZ
Washing Control	To Enter Water 2~8L/min Automatically
Products	Mobile Medical X Ray Machine
Kv Range	40-150 KV
Mas	0.4-360mAs
Rectifier Method	Full-wave Rectifier
Nominal Power	100kV, 40mA, 0.1s, 4kW

This portable X Ray machine is a perfect veterinary radiology equipment solution that is most suitable for small and medium-sized animal hospitals, pet hospitals, clinics, and other pet medical institutions. It can also be used as a digital x ray system for human dental imaging purposes.

Applications:

The VOLITION Medical X Ray Machine is a state-of-the-art digital radiography equipment that is suitable for a wide range of medical applications. This high-quality X-ray machine is designed and made in China with the highest quality standards and has been certified by ISO13485 and FSC.

The VOLITION Medical X Ray Machine is designed to operate on a power supply of 220V/50Hz, which makes it suitable for use in a variety of settings. It is equipped with a wide range of features and settings that make it ideal for use in different medical applications.

The VOLITION Medical X Ray Machine has a Mas range of 0.4-360mAs, which makes it suitable for use in a variety of imaging procedures. Additionally, it has an intermittent mode of 100kV, 200mA, and 20kW, which makes it ideal for use in industrial radiography equipment.

The VOLITION Medical X Ray Machine is equipped with a rotating anode X-ray tube that provides high-quality images with a high level of clarity and detail. This makes it suitable for use in a wide range of medical applications, including dental imaging, where the handheld dental X-ray machine is used.

The VOLITION Medical X Ray Machine has a maximum current of 500mA, which makes it ideal for use in a wide range of medical applications that require high levels of accuracy and precision. With its advanced features and reliable performance, the VOLITION Medical X Ray Machine is the perfect choice for medical professionals who require top-quality imaging equipment for their patients.

Overall, the VOLITION Medical X Ray Machine is the perfect choice for a range of medical applications, including dental imaging, industrial radiography equipment, and other medical procedures that require high-quality imaging equipment that provides accurate and precise results.

Customization:

Our VOLITION Mobile Medical X Ray Machine is a state-of-the-art digital x-ray system that is perfect for medical facilities. Our product is made in China and certified with ISO13485 and FSC. It uses a full-wave rectifier method and has a maximum current of 500mA. The voltage frequency is 380V/220V 50HZ and it has automatic washing control that enters water 2~8L/min.

Our customization services include:

Product branding with your company logo and name

Customized voltage frequency to meet your facility's requirements

Installation and training services for your staff

Customized washing control options for your convenience

Our product is not only limited to medical facilities, but can also be used for veterinary radiology equipment. Contact us today to learn more about our customization services and how we can help meet your needs.

Packing and Shipping:

Product Packaging:

The Medical X Ray Machine will be securely packed in a sturdy cardboard box.

The box will be labeled with the product name, model number, and handling instructions.

Inside the box, the machine will be wrapped in protective foam padding to prevent any damage during shipping.

Shipping:

The Medical X Ray Machine will be shipped via a reputable courier service.

Shipping costs will be calculated based on the destination and weight of the package.

The estimated delivery time will be provided to the customer upon completion of the order.

The customer will receive a tracking number to monitor the delivery status of the package.



Liuyang Volition Medical Equipment Co., Ltd.



+86 17775818268



258031353@qq.com



xray-system.com

76 Jianshou Ave, Economic and Technological Development Zone, Liuyang, Hunan, China