

OEM Veterinary Digital Radiography System Mobile X Ray Unit Radiology

Our Product Introduction

for more products please visit us on xray-system.com

Basic Information

- Place of Origin: China
- Brand Name: VOLITION
- Certification: ISO13485、FSC
- Minimum Order Quantity: 1
- Price: Negotiable
- Packaging Details: Standard Exportation Package
- Delivery Time: 7-15 work days
- Payment Terms: LC, T/T, D/P, PayPal, Western Union



Product Specification

- Exposure Time Range: 0.001-10 Seconds
- Image Resolution: ≥ 3.5 Lp/mm
- Imaging Technology: Digital
- Power Source: Battery
- Type: Portable
- Maximum Output Power: 100kV
- Maximum Output Current: 10mA
- Weight: Less Than 20kg
- Highlight: **OEM veterinary digital radiography system,
Mobile veterinary digital radiography system,
Radiography mobile x ray unit radiology**

Product Description

Product Description:

With a maximum output power of 100kV, this veterinary xray equipment is capable of producing high-quality x-ray images with excellent clarity and detail. The image resolution of this device is also impressive, with a minimum of ≥ 3.5 Lp/mm. This makes it ideal for diagnosing dental issues in animals of all sizes, from small pets to large livestock.

One of the standout features of this product is its power source. It is equipped with a rechargeable battery, which means that it can be used anywhere, even in remote locations, without the need for a power outlet. This makes it an excellent choice for veterinarians who need reliable and portable equipment for fieldwork or on-the-go diagnostic needs.

Whether you are a small animal veterinarian or work with large animals, this mobile digital x-ray machine is an essential tool for your practice. With its exceptional image quality, quick exposure time range, high output power, excellent image resolution, and reliable battery power source, it is the perfect choice for any veterinary xray equipment needs.

Features:

Product Name: Mobile Digital X Ray Machine

Maximum Output Current: 10mA

Image Detector: Flat Panel

Maximum Output Power: 100kV

Imaging Technology: Digital

Display Screen: LCD

This Mobile Digital X Ray Machine is perfect for a variety of applications, including use as a c arm fluoroscopy machine or veterinary xray equipment. It also works well as veterinary dental radiography equipment.

Technical Parameters:

Imaging Technology	Digital
Display Screen	LCD
Weight	Less Than 20kg
Certifications	CE, FDA
Maximum Output Power	100kV
Maximum Output Current	10mA
Type	Portable
Power Source	Battery
Exposure Time Range	0.001-10 Seconds
Image Detector	Flat Panel

This mobile fluoroscopy machine is a portable DR X Ray, suitable for dental X Ray devices.

Applications:

One of the primary applications of VOLITION mobile digital X-ray machines is dental radiographic equipment. These machines are designed to produce high-quality X-ray images of teeth, bones, and soft tissues in the mouth. They are also small and portable, making them easy to transport between dental offices and clinics. With exposure time ranges of 0.001-10 seconds, VOLITION mobile digital X-ray machines can produce high-quality images in a matter of seconds.

Another application of VOLITION mobile digital X-ray machines is in mobile fluoroscopy machines. These machines are used to produce real-time images of internal organs and tissues. They are commonly used during surgeries to guide the surgeon and ensure that the procedure is successful. VOLITION mobile digital X-ray machines are ideal for this application due to their real-time imaging capabilities and advanced imaging technology.

VOLITION mobile digital X-ray machines are also suitable for use in veterinary clinics. They can be used to produce high-quality X-ray images of animals, which is essential for diagnosing illnesses and injuries. These machines are small and portable, which makes them easy to transport to different locations.

With a maximum output power of 100kV and a flat panel image detector, VOLITION mobile digital X-ray machines produce high-quality images that are clear and accurate. They are also certified by CE and FDA, which ensures that they meet the highest standards for safety and quality.

Overall, VOLITION mobile digital X-ray machines are suitable for a wide range of applications in the medical field. They are ideal for dental radiographic equipment, mobile fluoroscopy machines, and veterinary clinics. With their advanced imaging technology and real-time imaging capabilities, VOLITION mobile digital X-ray machines are a valuable tool for medical professionals.

Customization:

Customize your VOLITION Mobile Digital X Ray Machine to fit your specific needs with our product customization services. Our dental radiographic equipment can be tailored to your requirements, whether you need a dental x ray machine panoramic or veterinary xray

equipment.

Our machines are manufactured in China and certified with ISO13485 and FSC, ensuring the highest quality. The machine operates on a battery power source, providing flexibility and mobility. The exposure time range is adjustable from 0.001 to 10 seconds, making it suitable for a wide range of applications. The machine is equipped with a flat panel image detector and can output up to 10mA maximum current. At less than 20kg in weight, it is lightweight and easy to transport.

Packing and Shipping:

Product Packaging for Mobile Digital X Ray Machine:

The Mobile Digital X Ray Machine will be shipped in a sturdy cardboard box that has been specifically designed to protect the machine during transportation. The box will be secured with packing tape and labeled with the necessary shipping information.

Inside the box, the Mobile Digital X Ray Machine will be packed securely with foam padding and other protective materials to ensure that it arrives at its destination in perfect condition.

Shipping:

The Mobile Digital X Ray Machine will be shipped via a reputable shipping company that offers fast and reliable delivery services. The shipping cost will be determined based on the destination and the preferred delivery method.

Customers will be provided with a tracking number so they can monitor the progress of their shipment. Delivery times may vary depending on the location, but we will do our best to ensure that the product is delivered as quickly as possible.



Liuyang Volition Medical Equipment Co., Ltd.



+86 17775818268



258031353@qq.com



xray-system.com

76 Jianshou Ave, Economic and Technological Development Zone, Liuyang, Hunan, China